

## Bechtol S plus ready-to-use solution

### Description

The highly effective disinfecting and cleaning solution for rotating instruments

### Application

Immediately after use place rotating instruments into Bechtol S plus. Heavily soiled instruments should be soaked in warm water then cleaned with a soft brush. Ensure that the instruments are totally covered by the solution. Avoid the formation of air pockets. Respect times of immersion. Then remove from the solution, rinse thoroughly and dry with a clean cloth. Observe instructions of manufacturer of instruments.

### Properties

- effective against bacteria (incl. MRSA) and yeast in 15 minutes (EN 13727 / EN 13624)
- guarantees active protection against tuberculosis (*M.terrae*) in 30 minutes (EN 14348)
- inactivates enveloped viruses (HBV/HCV/HIV/BVDV/Vaccinia) and non-enveloped viruses (Adeno/Polio/HPV/Rota and Noro) in 15 minutes (EN 14476)
- effective against spores in 30 minutes (EN 13704)
- contains no phenols or aldehydes
- suitable for use in an ultrasound device in 5 minutes
- excellent compatibility on materials



### Composition

100 g Bechtol S plus contain:

- 0,24 g Benzalkoniumchloride
- 0,12 g Dodecylpropylentriamine

### Shelf life

30 months

### Form of delivery

**REF 538** 1L bottle, ready-to-use solution

**REF 538.5** 5L canister, ready-to-use solution

### Classification

Medical device, class IIb, rule 15,  
CE Certificate: G1 17 09 49628 010,  
TÜV SÜD Product Service GmbH, München

### Manufacturer

Alfred Becht GmbH  
Carl-Zeiss-Str. 16  
D-77656 Offenburg  
Tel.: +49 781 60586-0  
Fax.: +49 781 60586-40  
Germany  
[www.becht-online.de](http://www.becht-online.de)

### Accessories

**REF 567.5** Dispenser pump for  
5 Liter canister

**REF 567.51** Pouring spot for 5  
and 10-Liter canister

**REF 630** Fräsator

**REF 631** Mini-Fräsator

