

Bechtozid soft

Ready-to-use Solution

Description

Foam-spray for the rapid disinfection of all alcohol sensitive medical devices, such as surfaces, treatment units, natural and synthetic leather, acrylics. For optimum hygiene and optimal patient flow.

Application

Spray Bechtozid soft on the contaminated surface. Respect time of immersion, then wipe surface with a dry and clean cloth. Wear suitable protective gloves and eye/face protection. The spray nozzle is adjustable to easily switch from a spray to a foam application.

Properties

- effective against bacteria (incl. MRSA) and yeast in 1 minute (EN 13624/13727)
- inactivates HBV/HCV/HIV/BVDV/Vaccinia-viruses in 1 minute (EN 14476)
- tested in accordance with the directives of the DGHM/VAH (Deutsche Gesellschaft für Hygiene und Mikrobiologie)
- contains no aldehydes or phenols
- biodegradable
- alcohol-free and without aerosol
- free of perfumes and coloring agents for best compatibility
- the canister is UN-approved and recyclable
- pH-value: 6,5



Composition

100 g of solution contain

- 250 mg Benzyl-C12-16-alkyldimethyl-Chloride
- 250 mg Didecyldimethylammonium-Chloride
- 250 mg C12-14 Alkyl[(ethylphenyl) methyl]dimethyl-Chloride

Form of delivery

- REF 537** 1-Liter bottle, ready-to-use solution, reusable, duo spray head
- REF 537.5** 5-Liter canister, ready-to-use solution
- REF 537.10** 10-Liter canister, ready-to-use solution

Classification

Medical device, class IIa, rule 15,
CE Certificate: G1 17 09 49628 010,
TÜV SÜD Product Service GmbH, München

Shelf life

30 months

Accessories

REF 537.SP Duo spray head for 1-Liter bottle

REF 567.5 Metering pump for 5 and 10-Liter canister

REF 567.51 Pouring spout for 5 and 10-Liter canister



Manufacturer

Alfred Becht GmbH
Carl-Zeiss-Str. 16
D-77656 Offenburg
Tel.: +49 781 60586-0
Fax.: +49 781 60586-40
Germany
www.becht-online.de