

Bechtifix alkoholfrei

Ready-to-use disinfecting wipes

Description

Recommended for disinfecting and cleaning of alcohol sensitive medical devices, natural and synthetic leather. Free of alcohol.

Application

Take a single tissue from the middle of the roll and wipe the surface thoroughly. Respect time of exposure. Then wipe surface with a dry and clean cloth. Close the box immediately after taking the tissue.

Properties

- effective against bacteria (incl. MRSA) and yeast in 1 minute under high organic load (EN 13727/13624)
- **successfully tested in the 4-field test acc. to EN 16615**
- inactivates HBV/HCV/HIV/BVDV/Vaccinia-viruses in 1 minute (EN 14476)
- particularly soft, absorbent and tear-resistant non-woven fabric
- free of perfumes and coloring agents for best compatibility
- high volume of release of disinfecting solution on surface
- disinfecting liquid: Bechtozid soft
- size of tissue: 14 x 20 cm



Composition

100 g of solution contains:

- 250 mg Benzyl-C12-16-alkyldimethyl-Chloride
- 250 mg Didecyldimethylammonium-Chloride
- 250 mg C12-14 Alkyl[(ethylphenyl) methyl]dimethyl-Chloride

Form of delivery

REF 549.100 Dispenser box with 100 tissues
REF 549.100NF Refill with 100 tissues

Expert opinions

Expertise according to DGHM:

- Prof. Dr. med. R. Schubert / Frankfurt am Main / 2010
- Prof. Dr. med Egon Marth / Graz / 2010

BVDV-/Vacciniavirus:

- Prof. Dr. Jochen Steinmann / Bremen / 2005

4-field-test (EN 16615):

- Prof. Dr. med. H.-P. Werner / Schwerin / 2015

Shelf life

30 months

Classification

Medical device, class IIa, rule 15,
CE Certificate: G1 17 09 49628 010,
TÜV SÜD Product Service GmbH, München

Manufacturer

Alfred Becht GmbH
Carl-Zeiss-Str. 16
D-77656 Offenburg
Tel.: +49 781 60586-0
Fax.: +49 781 60586-40
Germany
www.becht-online.de